UV 3100

4120021

**ESPECTROFOTÓMETRO
SPECTROPHOTOMETER**

Indice

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Index

1 Introduction	31
1.1 Packing List	31
1.2 Accessories	31
1.3 Appearance	31
1.4 Operation Panel	32
1.5 Keys Description	32
2 Safety	33
2.1 Safety icons	33
2.2 Risks to which the operator is exposed	34
2.3 Staff qualification	34
2.4 Modifications	34
3 Operating Principle	35
3.1 Main Functions	36
4 Installation	37
4.1 Location	37
4.2 Power supply connection	37
4.3 Installation	38
5 Operation	39
5.1 Basic Operation	40
5.2 Before Measurement	41
5.3 Measurement	41
5.4 System Utility	48
6. Equipment Maintenance	52
6.1 Daily Maintain	52
6.2 Troubleshooting	52
6.3 Spare parts replacement	53
7 Guarantee	56
8 Technical Specifications	56

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English Version

1 Introduction

UV-3100 Spectrophotometer has the features of wide range wavelength, high sensitivity, powerful functions, easy to use, simple and nice structure. Besides these, its large LCD, high precise A/D and easy to store RAM makes this equipment much more superior than other originals. It is widely used in chemistry, pharmaceuticals, biochemical, metallurgy, light industry, textile, medical, education, etc. It is one of the most important equipments in quality control and it is essential in normal laboratories.

Handle the parcel with care. Carefully unpack and check that the contents coincide with the packing list. If you see any part damaged or missing, please advise the distributor or J.P. Selecta, s.a. immediately.

1.1 Packing List

- Spectrophotometer **UV-3100**
- Power connection cable
- Instructions manual
- Communication cable
- Software management
- Equipment dust protection cover
- Glass cuvettes, a set of 4 units
- Quartz cuvettes, a set of 2 units

1.2 Accessories


- Thermostated flow Peltier cuvette holder 4120018
- Cuvette holder (path length 1-50mm) 4120030
(path length 100mm) 4120031
- Test tubes holder 4120032
- Ink printer (not thermal) 4120117

1.3 Appearance

1. LCD display
2. Keypad
3. Lid
4. Rod
5. USB port (memory)
6. LCD adjustment
7. Printer port
8. USB port (communication)
9. Cover of fan
10. Power socket
11. Power switch
12. Cover of cooling vents

1.4 Operation Panel


1.5 Keys Description

Function Keys, choose corresponding commands listed in the bottom of the LCD.

Note:

To describe conveniently, we call them **F1, F2, F3, F4** separately from left to right in this manual.

Input Parameters, Wavelength and set 8 cell auto cell holder

Delete the input value or stored data

Cancel operation or return to former me

Calibrate 100%T/0A

Load the files saved in the R

Save the test result

Start/Cancel Te

Set Waveleng

Print Test Resu

Confirm the input value or setti

Select/Release auto cell hold

Scroll the menu or da

Seek wavelength peak or set coordinate.

2 Safety

Spectrophotometers **UV-3100** include the safety measures suitable for a safety operation in a laboratory.

This equipment has been designed to be used by technicians who have been qualified for its use. For the correct and safety use of this equipment it is important that the laboratory staff know the safety general procedures, together with the safety precautions explained in this manual.

Some of the chemical products used in this spectrophotometer are corrosive and/or inflammable and the samples can be radioactive, toxic or potentially infectious. Be careful to follow the safety standardized procedures in laboratories for samples and chemical products manipulation.

The risky situations are indicated along this manual by means of the safety icons and together with some instructions that must be respected.

2.1 Safety icons

They identify risky situations and the safety measures to be taken.

The icons make reference to the paragraph marked in grey.

Danger risk

Risk of danger.

Please, respect the instructions indicated to do the described operation.

Electric risk, high voltage!

Risk of electric accident when accessing the areas labelled with this signal or when doing the operations indicated with this icon in the present manual.


Please, respect the instructions indicated to do the described operation.

Burn risk due to contact with high temperature areas

The temperature in the area indicated by this icon can exceed 60°C. Please, use thermal resistant gloves to do the described operation.

Please, respect the instructions indicated to do the described operation.

Ground

Fuse

Waste

2.2 Risks to which the operator is exposed

The working process with the described equipment exposes the operator to the following risky situations:

- Corrosive and / or inflammable substances manipulation.
- Sample manipulation which can be radioactive, toxic or potentially infectious.
- Glass pieces manipulation.
- Electric risk.

2.3 Staff qualification

This equipment can only be used by staff qualified in knowing dangers exposed in chemical analysis laboratories.

This equipment can only be used by conveniently trained staff who have read the instruction manual or have been qualified in the use of this equipment.

2.4 Modifications


The modification of the equipment safety systems operation or manipulation by personnel not authorized by the manufacturer can expose the operator to risks not foreseen in this manual.

3 Operating Principle

UV-3100 equipment is composed of five parts:

- An halogen lamp and a deuterium lamp as light sources.
- A monochromator to select the wavelength and to eliminate the unwanted second order radiations.
- A compartment to accommodate the sample cuvettes.
- A photodetector to receive light transmittance and to convert it to electrical sign.
- A display where absorbance, transmittance and concentration values are shown.

The halogen lamp light (or the deuterium one if we work in a wavelength lower than 340 nm) is focused on the monochromator entrance where a collimated mirror points the light beam towards the diffraction grating.

The diffraction grating diffuses the light beam to produce a spectrum. A part of it focuses in the monochromator output slot by means of another collimated mirror.

The wavelength light beam selected passes through a filter which helps eliminating unwanted second order radiations from the diffraction grating.

This light beam goes through the sample, and then the photodetector produces an electrical sign used to make the transmittance or absorbance calculations, shown in the digital display.

3.1 Main Functions

UV-3100 has 7 functions:

- **Photometry**

This is the basic mode, the photometric mode can be switched between Absorbance and Transmittance. You can use 2-points method to measure the concentration.

- **Quantitative**

Two methods are included in the menu: Coefficient Method and Standard Curve Method. The Regression Equation and the test result can be stored in the RAM and be printed out.

- **Wavelength Scan**

Scan range and scan intervals can be set (0.1nm, 0.2nm, 0.5nm, 1.0nm, 2.0nm, 5.0nm). The data and the curve can be stored and printed out.

- **Kinetics**

The minimum sampling interval is 0.5 second. The data and the curve can be stored and printed out.

- **DNA/Protein Test**

Provide two methods. Customers can also set the parameters. The data can be printed out.

- **Multi-wavelength Test**

Users can set several wavelengths before test, then the system will test all the T or A values of the setting wavelength one by one. The results can be printed out.

- **Utility**

Users can set the favourable items to keep the instrument work under best conditions.




Figure 2

4 Installation

Once the equipment is unpacked, check it is all right. Preserve the packing for some days.

4.1 Location

To ensure the best performance, the following conditions are required:

Place the equipment over a flat, levelled and stable surface, suitable for the equipment weight and not exposed to direct sunlight.

The best working temperature range is 16-35°C and the humidity is 45-80%.

With the aim of having the best equipment performance, keep it isolated, as much as possible, from any magnetic or electrical field source or any other electrical device that could generate high frequency fields.

Install the unit in an area free of dust, corrosive gas and hard vibrations.

Place the equipment at a distance lower than a 1.5m from a power socket.

To work in a comfortable way, left a free space of 50cm around the equipment and a minimum of 20cm behind it.

Remove obstructions or materials that could obstruct the air flow around and underneath the equipment.

4.2 Power supply connection

Please, select a socket near the equipment and suitable for its power.

Use only the supplied cord or similar ones to connect the equipment.

For your safety, the equipment must be earthed.

Before switching on, make sure that the voltage and the power frequency correspond to the indicated in the equipment.

If the local voltage is not stable enough, a manostat is required.

Be sure to set the instrument's voltage switch at your local power supply, or severe damage may occur!

4.3 Installation

Step 1: Check the packing list

Unpack the contents, check the materials with the Packing List. Any damage or something missing, please contact the local dealer or J.P.Selecta, s.a.

Step 2: Position

Place the instrument on a stable surface carefully.

Step 3: Install printer (optional)

Make sure the printer's power is off. Link the printer's data cable to the equipment's parallel port.

Step 4: Connect the power cord

Make sure the instrument's power switch is in the Off condition, connect the power cord to the equipment. And insert another end in the socket provided with a protective earth contact.

Step 5: Switch on the power

Check again. Make sure that all the cables are right. Switch on the power. Then the instrument begins to self-test. After the self-test and 15 minutes pre-warm, it can work. The self-test includes the following steps:

D2 Lamp Pre-warm --> Com-port Initialization --> Initialize Printer --> Start kernel --> Initialize AD --> Filter position --> Pre-warm.


After pre-warming, the instrument will ask the user to re-calibrate the system. Users can decide if they need to re-calibrate the system or not. After this step, the equipment will work normally.

After transportation or a period of use, the working condition of the equipment may changed because of environment changing or any other cause (attenuation of lights' energy, vibration, etc.) and then users need to re-calibrate the system to keep it in the best working conditions.

5 Operation

This is the software structure of **UV-3100**.

5.1 Basic Operation

5.1.1 Select test method

In the main menu, press numeric keys or use  or  keys to choose the corresponding menu, and then press **[ENTER]** to confirm.

5.1.2 Set wavelength

Press **[GOTO]** to set wavelength, use numeric keys to input the value and press **[ENTER]** to confirm. When the wavelength changes to the point you set, the system calibrates 100%T/0Abs automatically.

5.1.3 Set parameters

In the different menus, press the corresponding F (1-4) to go into different parameter setting screen, use  &  to choose or input the values by means of the numeric keys. Then press **[ENTER]** to confirm, and press **[ESC]** to return.

5.1.4 Set auto-cell holder (optional accessory)

Press **[CELL]** to active the auto cell holder and press the numeric keys (1-8) to make corresponding cell position at the light path. Press **[CELL]** again to inactive the auto cell holder.

5.1.5 Delete the input value

Press  to delete the last character, and press **[CLEAR]** to delete all the characters.

5.1.6 Delete the test result and stored data

In the test screen, press **[CLEAR]** to delete the test result or stored data.

5.1.7 Calibrate 100%T/0Abs

Put the reference in the light path, then press **[100%t/0 ABS]** to calibrate 100%T/0Abs.

5.1.8 Measure samples

In the test screen, put the samples to be tested in the light path, and press **[START/STOP]** to measure.

5.1.9 Print the test result

In the test screen, press **[PRINT]** to print the test result.

5.1.10 Store the test result

In the test screen, press **[SAVE]** to store the test result. Input the file name by means of the numeric keys and press **[ENTER]** to confirm.

5.1.11 Load the stored file

In the test screen, press **[OPEN]** to go into file selecting screen. Press  &  to choose the file you want, and press **[ENTER]** to open.

In different menus, users can only load the files with the corresponding extension. These are some of the corresponding extensions:

- Quantitative *.qua
- Standard Curve *.fit
- Wavelength Scan *.wav
- Kinetics *.kin
- DNA/Protein test *.dna
- Multi-wavelength *.mul

5.2 Before Measurement

5.2.1 Self-test

Remove all the blocks in the light path and close the lid of the compartment. Switch on the power supply to begin the self-test.

5.2.2 Warm up

After self-test, the instrument goes into pre-warm state. For accurate test, at least 30 minutes of warm up is required.

5.2.3 Check the Cuvettes

The cuvettes must be clear and there must be no remains of the samples on the surface of it. Only Silicon (Quartz) cuvettes are permitted to be used in the range of UV area.

5.3 Measurement

5.3.1 Basic Mode

- **Step 1, Go into the Test screen in Basic Mode:**

In the main menu, press the key of [1] or press \triangleup & \triangledown to choose «Basic Mode» followed with [ENTER] pressed to go into the test screen (Fig. 4.1).

- **Step 2, Set Photometric Mode:**

Press F2 to set photometric mode. Press \triangleup & \triangledown to choose «Abs.», «T%» or «Conc./Factor» and press [ENTER] to confirm.

If users choose «Abs.» or «T%», please go to step 5 directly.

- **Step 3, Set Concentration Unit:**

Press F1 to set concentration unit. Press \triangleup & \triangledown to choose unit followed with [ENTER] pressed to confirm. You can also choose «other» to input the self defined unit.

- **Step 4, Set «Factor» or «Standard»**

Two methods are under your choice:

Method 1: Input Factor F

Press F3 to set F. Input the value of F by numeric keypad, press [ENTER] to confirm. Then the F value would display on the LCD.

Method 2: Standards Mark

Put the Reference Sample in the light path and calibrate 100%T/0Abs. Put the standard sample in the light path, press F4 to start the mark. Input the concentration value of the standard and press [ENTER] to confirm, then it displays on the screen.

- **Step 5, Set Wavelength:**

Press [GOTO] to set wavelength, input the value by the numeric keypad followed with [ENTER] pressed to confirm.

- **Step 6, Calibrate 100%T/0Abs:**

Put the Reference in the light path and press [100%T/0 ABS] to calibrate 100%T/0Abs.

- **Step 7, Sample Measurement:**

Put the sample to be measured in the light path, then the result displays on the screen automatically.

- **Step 8, Print Data:**

Press [PRINT] to print the result.


Fig. 4.1


Fig. 4.5


Fig. 4.6


Fig. 4.7

5) Calibrate Standard Samples. Put the corresponding standard samples in the light path as the screen indicates and press **[START/STOP]** to measure. Then the Abs. value would appear in the corresponding table (Fig. 4.5)

Load the Stored Curves

In the «Calibration Table» screen, press **[OPEN]** to go into files select screen. Use & to select the curve you need and press **[ENTER]** to load.

Users can press F4 to see the curve (Fig. 4.6), and then press **[ESC]** to cancel.

• Step 4, Return the Sample Measurement Screen:

In the «Calibration Table» screen, press **[ESC]** to return to the Quantitative Test screen (Fig. 4.7).

• Step 5, Calibrate 100%T/0Abs:

Put the Reference Sample in the light path, press **[100%T/0 ABS]** to calibrate 100%T/0Abs.


Fig. 4.8

Step 6, Measure Samples:

Place the sample to be tested in the light path, press **[START/STOP]** to measure. Then the test result will be displayed in the data sheet. Repeat this step to finish measuring all the samples (Fig. 4.8).

Step 7, Print Data:

Press **[PRINT]** to print the result.

Step 8, Delete Data:

Press **[▲]** & **[▼]** to move the cursor of «*» to the value you don't need, then press **[CLEAR]** to delete the value.

Step 9, Store the Data:

After measurement, press **[SAVE]** to store the data. Input the file name by the numeric keypad and then press **[ENTER]** to store.


Fig. 4.9

5.3.3 Wavelength Scan

Step 1, Go into Wavelength Scan Screen:

In the main menu, press the numeric key **[3]** or **[▲]** & **[▼]** to choose «Wavelength Scan» and press **[ENTER]** to enter (Fig. 4.9).

Step 2, Parameters Setup:

Press F1 to set parameters, set «Scan From», «scan to», «scan step» and «scan speed», press **[ENTER]** to confirm.

Step 3, Set Photometric Mode:

Press F2 to set photometric mode, choose «T%», «Abs.» or «E» and press **[ENTER]** to confirm.

Step 4, Scan Baseline:

Place the Reference Sample in the light path, press **[100%T/0 ABS]** to scan the baseline, and press **[ESC]** to cancel.


Fig. 4.10

Step 5, Scan Samples:

Put the samples in the light path, press **[START/STOP]** to scan the sample (Fig. 4.10), press **[ESC]** to cancel.


Fig. 4.11

Step 6, Search Peaks:

After Scan, press F3 to go into peak search screen. Press F1 to set peak height, input the peak height and press [ENTER] to confirm. Press & to display the value of every wavelength point. Press & to display the value of every peak (Fig. 4.11).

Step 7, Smooth the Curve:

After scan, if there are many burrs on it, press F4 to smooth the curve.

Step 8, Print:

Press [PRINT] to print the curve and peaks value.

Step 9, Delete:

Press [CLEAR] to delete the local curve.

Step 10, Store the curve:

After scan, press [SAVE] to go into the store screen. Input the file name by means of the numeric keypad and then press [ENTER] to confirm.


Fig. 4.12

5.3.4 Kinetics

Step 1, Go into Kinetics Test screen:

In the main menu, press [4] or & to select «Kinetics Mode» and press [ENTER] to confirm (Fig. 4.12).

Step 2, Setup Parameters:

Press F1 to set parameters, input the corresponding values of «Total Time», «Delay Time» and «Time Intervals» according the screen indicates. Press [ENTER] to confirm.

Step 3, Set Photometric Mode:

Press F2 to set photometric mode, choose «T%» or «Abs.» and press [ENTER] to confirm.

Step 4, Set Wavelength:

Press [GOTO] to set wavelength, input the value of the wavelength by means of the numeric keypad and press [ENTER] to confirm.

Step 5, Calibrate 100%T/0Abs:

Put the Reference Sample in the light path, press [100%T/0 ABS] to calibrate 100%T/0Abs.

Step 6, Measure the sample:

Put the sample to be tested in the light path and press [START/STOP] to begin the test (Fig. 4-13), then press [ESC] to cancel.

Step 7, Calculate Response Rate:

After time scan, if users want to calculate the response rate of some period, you can


Fig. 4.13



Fig. 4.14

press F3 to go into «Process» screen. Input the values of «Begin Time», «End Time» and «Factor» separately and press **[ENTER]** to confirm. Then the value of «I.U.» would be displayed on the screen (Fig.4.14).

• **Step 8, Search:**

After scan finished, press F4 to go into search mode. Press  &  to search the value of every point (Fig. 4.15).

• **Step 9, Print:**

Press **[PRINT]** to print curve.

• **Step 10, Delete the Curve:**

Press **[CLEAR]** to delete the local curve.

• **Step 11, Save the Curve:**

After scan, press **[SAVE]** to save, input the file name by the numeric keypad. Press **[ENTER]** to confirm.



Fig. 4.15

5.3.5 DNA/Protein Mode

• **Step 1, Go into DNA/Protein Test:**

In the Main menu, press **[5]** or  &  to choose «DNA/Protein Mode» and press **[ENTER]** to confirm (Fig. 4.16).

• **Step 2, Setup Coefficients:**

Press F1 to set coefficients, input all the values of F1 to F4 by means of the numeric keypad according to the indication and press **[ENTER]** to confirm.

• **Step 3, Choose Measure Method:**

Press F2 to set method. Press  &  to choose «Absorbance Difference 1» or «Absorbance Difference 2» followed by **[ENTER]** pressed to confirm. If users don't want to measure reference, use  &  to choose «No» followed by **[ENTER]** pressed to confirm the choice.


Fig. 4.16

• **Step 4, Set Concentration Unit:**

Press F3 to set concentration unit. Use  &  to select unit and press **[ENTER]** to confirm.

• **Step 5, Calibrate 100%T/0Abs:**

Place the Reference Sample in the light path, press **[100%T/0 ABS]** to calibrate 100%T/0Abs.

• **Step 6, Sample Measure:**

Place the sample to be tested in the light path, press **[START/STOP]** to measure. The result will be displayed in the data sheet (Fig. 4-17).


Fig. 4.17


Fig. 4.18


Fig. 4.19

• **Step 7, Print Data:**
Press **[PRINT]** to print the result.

• **Step 8, Delete Data:**
Press **[CLEAR]** to delete the data.

• **Step 9, Store the Data:**
After measurement, press **[SAVE]** to store the result. Input the file name by the numeric keypad, then press **[ENTER]** to confirm.

5.3.6 Multi Wavelength Mode

• **Step 1, Go into Multi Wavelength Measurement:**
In the Main Menu, press **[6]** or **▲** & **▼** to choose «Multi Wavelength Measurement» and then press **[ENTER]** to go into multi wavelength measurement screen (Fig. 4.18).

• **Step 2, Setup Wavelength:**
Press F1 to go into wavelength setting screen, input all the wavelength value one by one by the numeric keypad. Press **[ENTER]** to confirm, press **[ESC]** to return.

• **Step 3, Set Photometric Mode:**
Press F2 to set the photometric mode, use **▲** & **▼** to select «Abs.» or «T%» mode, press **[ENTER]** to confirm.

• **Step 4, Calibrate 100%T/0Abs:**
Place the Reference Sample in the light path, press **[100%T/0 ABS]** to calibrate 100%T/0Abs.

• **Step 5, Measure Sample:**
Place the sample to be tested in the light path, press **[START/STOP]** to measure, the test result will be displayed in the data sheet (Fig. 4.19).

• **Step 6, Print Data:**
Press **[PRINT]** to print the result.

• **Step 7, Delete Data:**
Press **[CLEAR]** to delete the data.

• **Step 8, Store the Data:**
After measurement, press **[SAVE]** to store the result. Input the file name by the numeric keypad and press **[ENTER]** to confirm.

5.4 System Utility

In the main menu, press **[7]** or use **▲** & **▼** to select «System Utility» and press **[ENTER]** to confirm (Fig. 4.20)


Fig. 4.20

5.4.1 Wavelength Reset

After some time, the energy of lights would drop off and there will be some but slight influence on the test result. In this case, user's can reset the wavelength to compensate. We suggest the user's to reset the wavelength every one or two months.

Press [1] or use & to choose «Wavelength Reset», then press [ENTER] to begin the calibration. During the course, opening the lid of the compartment is prohibited.

5.4.2 Printer Setup

Press [2] or use & to select «Printer Setup» and then press [ENTER] to confirm (Fig. 4.21).

- **Reset Printer**

Press [1] or use , to choose «Reset Printer» and press [ENTER] to confirm. Then the printer will resume the initial condition.

- **Select print port**

Press [2] or use & to choose «Select print port» and press [ENTER] to confirm. Use & to choose «LPT» or «Comm» and then press [ENTER] to confirm.

- **Select Printer**

Press [3] or use & to choose «Select printer» and press [ENTER] to confirm. Use & to select the printer's model and then press [ENTER] to confirm.

- **Change Print Mode**

Press [4] or use & to choose «Change print mode» and press [ENTER] to confirm. Two modes are under your choice, one is «Print Data Sheet», and the other is «Print the display interface».


Fig. 4.21

5.4.3 Lamp Service

Press [3] or use & to choose «Lamp Service», press [ENTER] to go into «lamps management» screen (Fig.4.22).


Fig. 4.22

• **Switch On/Off D2 Lamp**

Users can switch Off the D2 Lamp when test in the wavelength range of 340-1100nm, so as to prolong the life time of D2 lamp.

Press [1] or & to choose «Switch On/Off D2 lamp» and then press [ENTER] to switch on or switch off the D2 lamp.

• **Reset the Usage Time of D2 Lamp**

After some period of usage (about 1000h), the energy would drop off to a low level, and the test data in the UV range would get instable and deviate a lot from the experiential value. Then the D2 lamp should be replaced to a new one. After replacement, users should calibrate the wavelength again.

Press [2] or use & to choose «Reset D2 lamp usage time» and press [ENTER]. Then users can find the D2 lamp's usage time and the system will ask if you are sure to reset the usage time. Press & to choose «Yes» and then press [ENTER] to confirm. The system will record the usage time from zero.

• **Switch On/Off W Lamp**

Users can switch Off the W Lamp when test in the wavelength range of 190-339nm, so as to prolong the life time of W lamp.

Press [3] or use & to choose «Switch On/Off W lamp» and then press [ENTER] to switch on or switch off the W lamp.

• **Reset the W Lamp Usage Time**

After some period of usage (about 1500h), the energy would drop off to a low level, and the test data in the Visible range would get instable and deviate a lot from the experiential value. Then the W lamp should be replaced to a new one. After replacement, users should calibrate the wavelength again.

Press [4] or use & to choose «Reset W lamp usage time» and press [ENTER]. Then users can find the W lamp's usage time and the system will ask if you are sure to reset the usage time. Press & to choose «Yes» and then press [ENTER] to confirm. The system will record the usage time from zero.

• **Change Switch Point**

The instrument permits the users to set the lamps changing wavelength point. In the range of 325—375, users can choose the lamp freely.

Press [5] or use  &  to choose «Change Switch point» and press [ENTER] to confirm. Input the wavelength point value (325—375nm) and press [ENTER] to confirm.

5.4.4 Clock Setup

Press [4] or use  &  to choose «Clock Setup», then press [ENTER] to go into «Clock Management» screen (Fig. 4.23).

- **Set Time**

Press [1] or  &  to choose «Set Time» and press [ENTER] to confirm. Input the time (Hour, Minute, Second) by means of the numeric keypad, press [ENTER] to confirm and return automatically.

- **Set Date**

Press [2] or use  &  to choose «Set Date» and press [ENTER] to confirm. Input the Date (Year, Month, Day) by numeric keypad, press [ENTER] to confirm and return automatically.

- **Display Time**

Press [3] or use  &  to choose «Display Time» and press [ENTER] to confirm. Then the time will be displayed on the Right Top corner.

- **Display Date**

Press [4] or use  &  to choose «Display Date» and press [ENTER] to confirm. Then the Date will display on the Right Top corner.


Fig. 4.23

5.4.5 Refresh Dark Current

When the working environment changed, users should refresh Dark Current before test.

Press [5] or use  &  to choose «Refresh Dark Current» and press [ENTER] to confirm. Then the system begins to refresh Dark Current.

Note: During the course, open the lid of the compartment is prohibited.

5.4.6 Connect to PC

Press [6] or use  &  to choose «Connect to PC», and press [ENTER] to go into waiting screen (Fig. 4.24). When the instrument was connected to PC, it displays «Controlled by PC».

5.4.7 Beeper On/Off

Press [7] or use  &  to choose «Beeper on/off» and press [ENTER] to switch on or switch off the beeper.

5.4.8 Language Selection


Fig. 4.24

Press **[8]** or use & to choose «Language Selection» and press **[ENTER]** to go into language selection screen, press & to choose «Chinese» or «English» and then press **[ENTER]** to confirm. Then, the screen will be changed into the language you just select.

5.4.9 Refresh System Baseline

Press **[9]** or use & to choose «Refresh System Baseline» and press **[ENTER]** to confirm. Then the system will refresh the baseline.

During this course, opening the lid of the compartment is prohibited.

5.4.10 Deleted Entire Saved Files

Use & to choose «Deleted Entire Saved Files» and press **[ENTER]** to confirm. Then the system ask you «Delete entire files, are you sure?» use & to choose «Yes», then all the files in the RAM will be deleted.

5.4.11 Restore Default Settings

Use & to choose «Restore Default Settings» and press **[ENTER]** to confirm, then the system will restore the initialization.

5.4.12 Setup Slit (Only for variable bandwidth models)

In the main menu, press **[8]** or use & to choose «Setup Slit» and press **[ENTER]** to confirm. Press & to choose «0.5nm», «1.0nm», «2.0nm» or «4.0nm» and then press **[ENTER]** to confirm.

6. Equipment Maintenance

To keep the instrument work in good condition, constant maintain is needed.

6.1 Daily Maintain

6.1.1 Check the Compartment

After measurement, the cuvettes with sample solutions should be taken out of the compartment in time. Or the volatilization of the solution would make the mirror go moldy. Users must pay more attention to the corrosive sample and liquid easy to volatilize. Any solution remains in the compartment should be wiped off immediately.

6.1.2 Surface Clean

The cover of the instrument is with paint. Please use wet towel to wipe off the drips on the surface immediately. Organic solution is forbidden to be used to clean the cover. Please wipe off the dirt on the cover timely.

6.1.3 Clean the Cuvettes

After every test or after a solution change, the cuvettes should be cleaned carefully, or the remains on the surface would cause measuring error.

6.2 Troubleshooting

6.2.1 Dark Current Error when Self-test

Possible Cause	Solution
Open the compartment lid during the course of self-test.	Close the compartment lid and switch on the power again.

6.2.2 No Response After Power On

Possible Cause	Solution
Bad contact in power supply Fuse melt	Improve the contact Replace a new fuse

6.2.3 Printer Not Work, Printer Error

Possible Cause	Solution
No power supply supply	Switch on the power
Bad contact in power supply Bad contact of the data cable	Improve the contact Improve the contact

6.2.4 No Stable Reading

Possible Cause

No enough pre-warm
 Glass cuvettes used in UV Range
 No stable Sample
 Much higher sample concentration
 Low voltage or unstable power supply
 Lights defect
 Light used up

Solution

Increase the pre-warm time
 Use Silicon Cuvettes.
 Improve the sample
 Dilute the sample
 Improve the power condition
 Replace a new lamp
 Replace a new lamp

6.2.5 Worse Repeatability

Possible Reason

Unstable sample
 Cuvettes polluted

Solution

Improve the sample
 Clean the cuvettes

6.2.6 Incorrect Reading

Possible Reason

Dark Current Error
 Worse matching of the cuvettes

Solution

Re-get the dark current
 Improve the matching of the cuvettes

6.3 Spare parts replacement

6.3.1 Replace the Fuse

Danger! Be sure to switch off the power and unplug the socket before replacement!


- **Step 1, Tools Preparation:**
Prepare a 3×75 Flat Blade screwdriver.

- **Step 2, Switch Off the Power Supply:**
Switch off the power supply, and unplug the socket.

- **Step 3, Take Out the Fuse Seat:**
Take out the fuse seat by the screwdriver (Fig. 5.1).

Fig. 5.1

- **Step 4, Replace a New Fuse:**
Pick out the spare fuse and replace it to the working position (Fig. 5.2).

- **Step 5, Reset the Fuse Seat:**
Replace the fuse seat in the power socket.

- **Step 6, Switch On the Power:**
Plug the socket and switch on the power.

Fig. 5.2

6.3.2 Replace Lamps

High temperature!! Wait 20 minutes before open the lamp chamber after power off to avoid scald!

- **Step 1, Tools Preparation:**

Prepare a 6×150mm Cross Blade screwdriver and a pair of glove.

- **Step 2, Power Off:**


Switch off the power supply and unplug the socket.

- **Step 3, Open the Cover:**

Unscrew the 4 screws indicated in Fig. 5.3 (each side with 2 screws) and remove the cover.

Fig. 5.3

- **Step 4, Open the Cover of the Light Chamber:**

Unscrew the 2 screws on the light chamber cover and remove it (Fig. 5.4).

Fig. 5.4


Fig. 5.5

- **Step 5, Replace the D2 Lamp:**

Unscrew the 2 screws on the D2 Flange (No.1 in Fig. 5.5), unplug the connector in the power board (No. 2) and remove the D2 lamp. Draw on the cotton glove and replace a new lamp. Fix the 2 screws and plug the connector again.

- **Step 6, Replace W Lamp:**

Remember the direction of the filament before pull out the W lamp. Be sure that the new lamp's filament is in the same direction as before.

Pull out the defected W lamp and draw on the cotton glove. Insert the new W lamp as deep as possible on the lamp seat. Be sure to keep the filament in the same direction as the old one face.




Fig. 5.6

• **Step 7, Adjust the Position of the W Lamp:**

Switch on the power (the switch mirror should be placed to the position as Fig. 5-7 indicates). Observe the entrance light, and it should be in the centre of the entrance hole (Fig.5-7). If the light deviates to left or right, then loosen the No.1 screws in Fig. 5-8 and move the lamp seat to the left or to the right until it focusses on the centre of the slot. Then fix the screws. If the light deviates to up and down, then loosen the No.2 screws in Fig. 5-8 and move the lamp seat up and down until the light focusses on the centre of the slot. Then fix the No. 2 screws again.

• **Step 8, Finish:**

Reset the cover of the light chamber and fix the screws. Reset the cover of the instrument and fix the screws. Recover the pole in the compartment, and then the course is finished.

Fig. 5.8

6.3.3 Replace the Battery

Be sure to switch off the power supply and unplug the socket before opening the bottom cover

• **Step 1, Prepare the Tools:**

Prepare a 6x150mm Cross Blade Screwdriver.


• **Step 2, Switch Off the Power Supply:**

Switch off the power supply and unplug the socket.

• **Step 3, Open the Bottom Cover Plate:**

Unscrew the 13 screws indicated in Fig.5-9 and then remove the bottom plate.

Fig. 5.9

• **Step 4, Replace the Battery:**

Pick out the old battery and replace a new one (Fig. 5-10).

• **Step 5, Finish:**

Recover the bottom plate and fix the 13 screws, and then the course id finished.

Fig. 5.10

7 Guarantee

This product is guaranteed for one year from the date of purchase and it covers defective components and assembly.

Please keep the purchase ticket as it is the only valid document for the purchase date.

This guarantee does not cover damage caused by an incorrect use of the equipment. This equipment must be properly validated before its use in the value sample analysis for the user.

Any return without the authorization of J.P. SELECTA, s.a. is not accepted.

8 Technical Specifications

Optical System	Single Beam
Wavelength Range	190-1100nm
Stray Light	=<0.05%T @ 220nm & 360nm
Band Width	2nm (UV-3100)
Photometric Range	0-200%T, -0.3-3.0A
Wavelength Accuracy	±0.5nm (UV-3100), =<0.3nm (UV-3100)
WL. Repeatability	±0.3%T (UV-3100)
Photometric Accuracy	±0.3%T (UV-3100)
Stability	0.002A/h @ 500nm (UV-3100)
Display	320×240 Dots Matrix LCD
Data Output	USB, Parallel Port
Dimension	490×360×240 (UV-3100)
Weight	14kg (UV-3100)