Referencia: PA6061

Ficha Técnica

Producto: Legionella Selective (GVPC) Agar

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- ATLAS, R.M. & L.C. PARKS (1993) Handbook of Microbiological Media. CRC Press. BocaRaton. Fla. USA.
- CLESCERI, L.S., A.E. GREENBERG & A.D. EATON (1998) Standard methods for the examination of water and wastewater. 9-106. 20th edition. APHA-AWWA-WEF. Washington DF, USA.
- EDELSTEIN, P.H., (1981) Improved semiselective medium for the isolation of Legionella pneumoniae from contaminated clinical and environmental specimens. J. Clin Microbiol. 14(3):298.
- FEELEY, J.C., R.J. GIBSON, G.W. GORMAN, N.C. LANGFORD, J.K. RASHEED, C.D. MACKEL, & W.B. BAINE (1979) Charcoal-Yeast Extract Agar: Primary isolation medium for Legionella pneumophila. J. Clin. Microbiol. 10(4) 437.
- ISO 11731 Standard (1998) Water Quality - Detection and Enumeration of Legionella.
- ISO 11731-2 Standard (2004) Water Quality - Detection and Enumeration of Legionella - Part 2: Direct membrane filtration method for waters with low bacterial counts. Williams & Wilkins, Baltimore, Md, USA.
- ISO/TS 11133-1: 2009. Microbiology of food and animal feeding stuffs.- Guidelines on preparation and production of culture media. Part 1: General guidelines on quality assurance for the preparation of culture media in the laboratory.
- ISO/TS 11133-2: 2003 Corr. 2004. Microbiology of food and animal feeding stuffs.- Guidelines on preparation and production of culture media. Part 2: Practical guidelines on performance testing of culture media.
- MacFADDIN, J.F. (1985) Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria.
- PASCULLE, A.W., J.C. FEELEY, R.J. GIBSON, L.G. CORDES, R.L. MYEROWITZ, C.M. PATTON, G.W. GORMAN, C.L. CARMACK, J.W. EZZELL & J.N. DOWLING (1980) Pittsburgh pneumonia agent: Direct isolation from human lung tissue. J. Infect. Dis., 141:727.
- WARD, K.W. (1995) Processing and interpretation of specimens for Legionella spp. In "Clinical Microbiology Procedures Handbook" Chap. 12.1 edited b H.D. Isenberg. ASM Press. Washington DC, USA.