

GENERAL INFORMATION

Product name : Crystallizing dish, high volume, polypropylene

Description : Without spout. Autoclavable at 121 °C

TECHNICAL DATA

reference	capacity	Ø (mm)	h (mm)	pcs/pack
CDIP-1K5-001	1500 ml	180	60	1
CDIP-3K0-001	3000 ml	200	100	1
CDIP-4K9-001	4900 ml	250	100	1

PACKAGING AND LOGISTICS

reference	vol (l)	kg	TARIC	GTIN
CDIP-1K5-001	2,0535	0,115	39233010	08434868012967
CDIP-3K0-001	4	0,185	39233010	08434868012974
CDIP-4K9-001	6,76	0,285	39233010	08434868012981

PRODUCT PHOTO



MATERIAL

MATERIAL Polypropylene

Good chemical resistance and autoclavable at 121°C.

INFORMACIÓN GENERAL

Nombre del producto : Cristalizador de gran volumen de polipropileno**Descripción :** Sin pico. Autoclavable a 121 °C

DATOS TÉCNICOS

referencia	capacidad	Ø(mm)	h (mm)	unidades por ref.
CDIP-1K5-001	1500 ml	180	60	1
CDIP-3K0-001	3000 ml	200	100	1
CDIP-4K9-001	4900 ml	250	100	1

EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (l)	kg	TARIC	GTIN
CDIP-1K5-001	2,0535	0,115	39233010	08434868012967
CDIP-3K0-001	4	0,185	39233010	08434868012974
CDIP-4K9-001	6,76	0,285	39233010	08434868012981

FOTO DEL PRODUCTO



MATERIAL

MATERIAL Polipropileno

Buena resistencia química y autoclavable hasta 121°C.

INFORMATIONS GÉNÉRALES

Nom produit : Cristallisoir, grand volume, en polypropylène

Description : Sans bec. Autoclavable à 121 °C

DONNÉES TECHNIQUES

référence	capacité	Ø (mm)	h (mm)	unités par ref.
CDIP-1K5-001	1500 ml	180	60	1
CDIP-3K0-001	3000 ml	200	100	1
CDIP-4K9-001	4900 ml	250	100	1

EMBALLAGE ET LOGISTIQUE

référence	vol (l)	kg	TARIC	GTIN
CDIP-1K5-001	2,0535	0,115	39233010	08434868012967
CDIP-3K0-001	4	0,185	39233010	08434868012974
CDIP-4K9-001	6,76	0,285	39233010	08434868012981

PHOTO PRODUIT



MATÉRIEL

MATÉRIEL Polypropylène

Bonne résistance chimique et autoclavable à 121 °C.

INFORMAZIONE GENERALE

Nome del prodotto : Cristallizzatore in polipropilene per grandi volumi**Descrizione :** Senza becco. Autoclavabile a 121 °C

DATI TECNICI

referenza	capacità	Ø(mm)	h (mm)	unità per ref.
CDIP-1K5-001	1500 ml	180	60	1
CDIP-3K0-001	3000 ml	200	100	1
CDIP-4K9-001	4900 ml	250	100	1

IMBALLAGGIO E DATI LOGISTICI

referenza	vol (l)	kg	TARIC	GTIN
CDIP-1K5-001	2,0535	0,115	39233010	08434868012967
CDIP-3K0-001	4	0,185	39233010	08434868012974
CDIP-4K9-001	6,76	0,285	39233010	08434868012981

FOTO DEL PRODOTTO



MATERIALE

MATERIALE Polipropilene

Buona resistenza chimica e autoclavabile fino a 121 °C.