

GENERAL INFORMATION

Product name : Conical head clip, steel

Description : Traditional clips made of steel wire


TECHNICAL DATA

| reference | type | pcs/pack |
|--------------|------------------|----------|
| CHCS-014-010 | for 14/23 joints | 10 |
| CHCS-019-010 | for 19/26 joints | 10 |
| CHCS-024-010 | for 24/29 joints | 10 |
| CHCS-029-010 | for 29/32 joints | 10 |

PACKAGING AND LOGISTICS

| reference | vol (l) | kg | TARIC | GTIN |
|--------------|---------|-------|----------|----------------|
| CHCS-014-010 | 0,351 | 0,031 | 82032000 | 08434868013810 |
| CHCS-019-010 | 0,351 | 0,031 | 82032000 | 08434868013827 |
| CHCS-024-010 | 0,48 | 0,046 | 82032000 | 08434868013834 |
| CHCS-029-010 | 0,6 | 0,081 | 82032000 | 08434868013841 |

PRODUCT PHOTO



INFORMACIÓN GENERAL

Nombre del producto : Pinzas de acero para esmerilados**Descripción :** Pinzas tradicionales de alambre de acero


DATOS TÉCNICOS

| referencia | tipo | unidades por ref. |
|--------------|------------------|-------------------|
| CHCS-014-010 | para bocas 14/23 | 10 |
| CHCS-019-010 | para bocas 19/26 | 10 |
| CHCS-024-010 | para bocas 24/29 | 10 |
| CHCS-029-010 | para bocas 29/32 | 10 |

EMBALAJE Y DATOS LOGÍSTICOS

| referencia | vol (l) | kg | TARIC | GTIN |
|--------------|---------|-------|----------|----------------|
| CHCS-014-010 | 0,351 | 0,031 | 82032000 | 08434868013810 |
| CHCS-019-010 | 0,351 | 0,031 | 82032000 | 08434868013827 |
| CHCS-024-010 | 0,48 | 0,046 | 82032000 | 08434868013834 |
| CHCS-029-010 | 0,6 | 0,081 | 82032000 | 08434868013841 |

FOTO DEL PRODUCTO



INFORMATIONS GÉNÉRALES

Nom produit : Pince de rodage inox**Description** : Pince traditionnelle en fil d'acier


DONNÉES TECHNIQUES

| référence | type | unités par ref. |
|--------------|----------------|-----------------|
| CHCS-014-010 | pour col 14/23 | 10 |
| CHCS-019-010 | pour col 19/26 | 10 |
| CHCS-024-010 | pour col 24/29 | 10 |
| CHCS-029-010 | pour col 29/32 | 10 |

EMBALLAGE ET LOGISTIQUE

| référence | vol (l) | kg | TARIC | GTIN |
|--------------|---------|-------|----------|----------------|
| CHCS-014-010 | 0,351 | 0,031 | 82032000 | 08434868013810 |
| CHCS-019-010 | 0,351 | 0,031 | 82032000 | 08434868013827 |
| CHCS-024-010 | 0,48 | 0,046 | 82032000 | 08434868013834 |
| CHCS-029-010 | 0,6 | 0,081 | 82032000 | 08434868013841 |

PHOTO PRODUIT



INFORMAZIONE GENERALE

Nome del prodotto : Clips in acciaio per giunti conici smerigliati**Descrizione** : Clips tradizionali in filo di acciaio

DATI TECNICI

| referenza | tipo | unità per ref. |
|--------------|------------------|----------------|
| CHCS-014-010 | per conici 14/23 | 10 |
| CHCS-019-010 | per conici 19/26 | 10 |
| CHCS-024-010 | per conici 24/29 | 10 |
| CHCS-029-010 | per conici 29/32 | 10 |

IMBALLAGGIO E DATI LOGISTICI

| referenza | vol (l) | kg | TARIC | GTIN |
|--------------|---------|-------|----------|----------------|
| CHCS-014-010 | 0,351 | 0,031 | 82032000 | 08434868013810 |
| CHCS-019-010 | 0,351 | 0,031 | 82032000 | 08434868013827 |
| CHCS-024-010 | 0,48 | 0,046 | 82032000 | 08434868013834 |
| CHCS-029-010 | 0,6 | 0,081 | 82032000 | 08434868013841 |

FOTO DEL PRODOTTO

