


Specification

**Pepsin 1:10.000 NF**

**175208**


<b>Physical Description:</b>	Solid
<b>Product Code:</b>	175208
<b>Product Name:</b>	Pepsin 1:10.000 NF
<b>Specifications:</b>	Proteolytic Activity: 1 : 10000 Insoluble matter in H <sub>2</sub> O: passes test Identity (IR): passes test pH (5 %; H <sub>2</sub> O): 3 - 5
<b>Hazard pictograms</b>	
<b>WGK:</b>	1
<b>Storage:</b>	2 - 8°C protected from light
<b>Signal Word:</b>	Danger
<b>GHS Symbols:</b>	GHS07 GHS08
<b>H Phrases:</b>	H315 H319 H334 H335
<b>P Phrases:</b>	P284 P305+P351+P338 P342+P311 P403+P233 P405

**AppliChem GmbH**

Specification

**Pepsin 1:10.000 NF**

**175208**

	P501
<b>CAS:</b>	9001-75-6
<b>EINECS:</b>	232-629-3
<b>CS:</b>	35079090
<b>Index Nr.:</b>	647-008-00-6
	

**AppliChem GmbH**

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • [info.de@itwreagents.com](mailto:info.de@itwreagents.com) • [www.itwreagents.com](http://www.itwreagents.com)  
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340