

# Technical Data Sheet

GRADE:	F1095
DESCRIPTION:	General purpose qualitative filter paper
DATE:	October 2022

The quality control data given in this document represent the quality of the released lot. These values are the basis for the official release of this material. The Quality Department for quality control of filters has measured the values and assures that they are within the limits that are established in the current specification for this material. The values stated do not represent any internal or external specification for this particular material. This product has passed external-house tests and thus meets Chmlab Group stringent quality control standards. The following is checked on a regular basis:

## Technical Specifications

PACKAGING:	Boxes of 100 units
FORMAT:	Circles (Ø mm) 25, 40.5, 42.5, 47, 50, 55, 70, 90, 100, 110, 125, 130, 150, 185, 200, 240, 250, 270, 300, 320, 350, 400, 450, 500. Other sizes under request
Basis weight	90 g/m <sup>2</sup>
Thickness	330 µm
Retention range	26-34 µm
Filtration speed	Medium
Ash content	< 0.1 %