

# Technical Data Sheet

|              |  |
|--------------|--|
| GRADE:       | F1096                                    |
| DESCRIPTION: | General purpose qualitative filter paper |
| DATE:        | October 2022                             |

The quality control data given in this document represent the quality of the released lot. These values are the basis for the official release of this material. The Quality Department for quality control of filters has measured the values and assures that they are within the limits that are established in the current specification for this material. The values stated do not represent any internal or external specification for this particular material. This product has passed external-house tests and thus meets Chmlab Group stringent quality control standards. The following is checked on a regular basis:

## Technical Specifications

|                  |   |
|------------------|---|
| PACKAGING:       | Boxes of 100 units  |
| FORMAT:          | Circles (Ø mm) 25, 40.5, 42.5, 47, 50, 55, 70, 90, 100, 110, 125, 130, 150, 185, 200, 240, 250, 270, 300, 320, 350, 400, 450, 500.<br>Other sizes under request |
| Basis weight     | 85 g/m <sup>2</sup>   |
| Thickness        | 180 µm  |
| Retention range  | 14 - 22 µm  |
| Filtration speed | Medium - fast   |
| Ash content      | < 0.1 %   |